



Hardness (Total & Calcium) Test Kit, Model HA-4P

Ürün Numarası:

145700

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Hardness, Total & Calcium drop count titration/EDTA test kit, range 17.1 - 342 mg/L CaCO₃, 100 tests.

The kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field.

Utilises drop count titration - a simple, economical titration method

Teknik özellikler

Kutudakiler:	Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.
Metot:	Drop Count Titration/EDTA
Muhafaza Türü:	D
Ölçüm aralığı:	17,1 - 342 mg/L CaCO ₃
Parametre:	Hardness, Total & Calcium - As CaCO ₃
Platform:	Drop Count Titration
Smallest increment:	17.1 mg/L
Test kiti modeli:	HA-4P
Test sayısı:	100

Kutudakiler

Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.