



Test tube/ MPN rack, polyethylene, 20 mm (O.D.), 40 holes

Ürün Numarası:

2497902

TRY Fiyat:
Mevcut

Bize Ulaşın

Forty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Teknik özellikler

Çap:	20 mm
Hacim / Paket boyutu:	1 each
Test tüpü düzenlemesi:	4 x 10
Test tüpü yuvası:	40