



## Bottle rack, plastic-coated, 20 compartments

Ürün Numarası: 2551801  
TRY Fiyat: Bize Ulaşın  
Mevcut

Plastic-coated wire rack provides corrosion resistant transport system for sample bottles and bags. Racks feature 15 or 20 compartments with double-grid bottom to prevent loss of samples during transport. 15 compartment rack is ideal for Whirl-Pak sample bags (4 oz). 20 compartment rack fits most bottles associated with standard drinking water sampling. Autoclavable.

### Teknik özellikler

Hacim / Paket boyutu: 1  
Önemli özellik: Rack, 20 Bottle Compartments