



**Be Right™**



## Test strips pH , 4-9 pH units

Ürün Numarası:

2745650

TRY Fiyat:

Bize Ulaşın

Mevcut

pH test strips offer a convenient and versatile solution for measuring the acidity or alkalinity of various aqueous solutions. Most natural waters exhibit pH levels ranging from pH 4 to pH 9, with common values exceeding pH 7 due to the presence of carbonates and bicarbonates, contributing to alkalinity. In industrial settings, maintaining precise pH levels in processes such as manufacturing, chemical production, and wastewater treatment is crucial for optimal efficiency and product quality. Similarly, in boiler feed water systems, controlling pH within a narrow range is essential to prevent corrosion and scale buildup. Additionally, pH regulation in swimming pools is vital for ensuring water comfort, swimmer safety, and the effectiveness of sanitizers.

With 50 tests per package, our pH test strips provide a cost-effective and user-friendly solution for accurate pH measurement in a wide range of applications.

**Easy to use, disposable, and inexpensive**

**Great for quick evaluations of pH levels in the lab or in the field**

### Teknik özellikler

Birim:	50 tests
En küçük artış adımları:	Steps: 4, 5, 6, 7, 8, 9 pH
Ölçüm aralığı:	4 - 9 pH units
Parametre:	pH, mid range
Platform:	Test Strip
Test kiti modeli:	Test Strips