



Digital titrator assembly

Ürün Numarası: 1690002
TRY Fiyat: Bize Ulaşın
Mevcut

Portable and robust digital Titrator for fast and precise titrations. High resolution (125µL/digit) and high accuracy ($\pm 1\%$). Easy handling. Applications for more than 40 parameters with interchangeable cartridges.

Accurate results - when and where you need them

Fast and Flexible

Versatile

Economical

Teknik özellikler

Aksesuarlar dahil: Hayır
İçerik: Digital Titrator device only