



Test Kit, chlorine (free & total) & pH, model CN-67

Ürün Numarası:

1411100

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Control of chlorine and pH in pools and spas to maintain safe bathing conditions and to protect against system corrosion. Supplied with color discs and comparator for monitoring chlorine and pH, the Free & Total Chlorine and pH Test Kit provides more accuracy than test kits using color cubes or step-wise comparators.

One of Hach's most popular test kits

Includes color discs and a comparator for monitoring chlorine and pH

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Model:	CN-67
Ölçüm aralığı:	6.6 - 8.4 pH units
Ölçüm aralığı (2):	0 - 3,4 mg/L
Parametre:	Chlorine (Free & Total) and pH
Platform:	Renk Diski
Test sayısı:	300
Uygulama:	Pool and Spa