



## Deluxe water conditioning demonstration test kit, model AD-16E

Ürün Numarası:

148410

TRY Fiyat:  
Mevcut

Bize Ulaşın

Tehlikeli



Designed to boost your testing capabilities, Hach's most complete water conditioning test kit includes tests and demonstrations for seven critical water quality parameters. A optional softener regeneration test kit is available (P/N - 2666000).

**Uses color disc and drop count titration methods**

**In addition to soap consumption and softener demonstrations, includes tests to measure iron, pH, and nitrate**

**Also includes the Total Dissolved Solids (TDS) Pocket Pal<sup>™</sup> Tester**

**Comes complete with a sturdy carrying case and features viewing tube, water softener, and a compartmentalized reagent cabinet to keep chemistries neat and well-organized**

---

### Teknik özellikler

Model:	AD-16E
Platform:	Renk Diski
Uygulama:	Demonstration