



Test kit, chlorine, hardness, iron & pH , model CN-39WR

Ürün Numarası:

223002

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Chlorine, Hardness, Iron, and pH Test Kit combines the most popular water treatment test into one compact kit. The chlorine test is USEPA-accepted for drinking water testing if approved by individual states.

Includes frequently measured water quality parameters in one rugged test kit

Uses simple, economical drop count titration and color disc methods

Tests chlorine (0-3.4 mg/L), hardness (1-20 gpg), iron (0-4 mg/L), and pH (4-10).

Teknik özellikler

Model:	CN-39WR
Ölçüm aralığı:	0 - 3,4 mg/L
Parametre:	Chlorine, Hardness, Iron, and pH
Platform:	Renk Diski
Test sayısı:	100
Uygulama:	Water Quality