



Test kit, professional boiler treatment/boilerfeed & cooling water, model PBC-DT

Ürün Numarası:

2350700

TRY Fiyat:

Bize Ulaşın

Teslimat tarihi için arayın

Tehlikeli



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. In addition to all of the tests offered in Model BT-DT (PN 2350600) this kit provides testing capability for iron, silica, chlorine, nitrite, and total dissolved solids (TDS). The versatility and accuracy of this kit will prove invaluable to boiler operators and chemical sales professionals.

Uses simple, economical digital titrator titration and color disc methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Platform:	Dijital Titratör
Test sayısı:	50
Uygulama:	Boiler Water