



Be Right™



Bottle, Wash, Nalgene, Narrow Mouth, 500 mL, Red Cap/ Stem, 6/pk

Ürün Numarası:

2928204

TRY Fiyat:
Mevcut

Bize Ulaşın

Narrow mouth polypropylene screw closure and stem feature one-piece construction to ensure a leak proof seal. Tip can be cut back to increase flow. 500 mL (16.8 oz. Capacity, 28 mm mouth). Choose from 4 cap/stem colors.

Teknik özellikler

Birim:	6 pk
Hacim / Paket boyutu:	6 pk
Kapasite:	500 mL
Malzeme:	LDPE
Renk:	Berrak
Stem:	Red Cap/Stem