



Cylinder, graduated, mixing, glass, 50 mL ± 0.4 mL, 1.0 mL divisions, glass stopper #16

Ürün Numarası: 189641
TRY Fiyat: Bize Ulaşın
Mevcut

Tall form, single metric scale, glass stopper, graduated. Calibrated "to contain," durable white scale.

Scale is durable white ceramic enamel. Designed to meet ASTM Specification E1272, Style 2 (Class B). Includes glass stopper.

Teknik özellikler

Gövde malzemesi:	Glass
Hacim / Paket boyutu:	1 each
Kapasite:	50 mL
Ölçüm metodu:	To Contain
Sertifikalar:	Class B
Tür:	Tek Metrik, Karıştırma