



Cylinder, graduated, mixing, Glass, 50 mL ± 0.4 mL, 1.0 mL divisions, glass stopper #16

Ürün Numarası:

2088641

TRY Fiyat:
Mevcut

Bize Ulaşın

This tall form glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes glass penny head stopper. Graduations are in durable white enamel bake.

Teknik özellikler

Gövde malzemesi:	Glass
Hacim / Paket boyutu:	1 each
Kapasite:	50 mL
Ölçüm metodu:	To Contain
Sertifikalar:	Class B
Tür:	Tek Metrik, Karıştırma