



## Bottle, wash, diamond alert, 500 mL, Variety 6/pk

Ürün Numarası:	2613049
TRY Fiyat:	Bize Ulaşın
Mevcut	

Wide-mouthed and color-coded wash bottles provide safe and positive identification of common laboratory solutions. Standard Diamond Alert symbols indicate type of hazard contents may represent. Polyethylene bottles with color-coded polypropylene screw caps and spray-tip assembly. Variety 6-pack includes one bottle each for water, acetone, methanol, isopropanol, toluene, and ethanol.

### Teknik özellikler

Birim:	6 pk
Hacim / Paket boyutu:	6 pk
Kapasite:	500 mL
Malzeme:	LDPE
Renk:	Transparent
Stem:	Blue Cap/Stem