



Test kit, iron & pH color disc , model IR-23

Ürün Numarası:

256200

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Utilizing the simple, accurate color disc method, this Iron & pH Water Quality Test Kit provides accurate, on-site measurement of two important parameters to help maintain water quality within the desired range. Range: 0-10 mg/L/ 4-10 pH unit.

Uses color disc to provide accurate, on-site measurement of two important parameters to help maintain water quality within the desired range

Teknik özellikler

Model:	IR-23
Parametre:	Iron and pH
Platform:	Renk Diski
Uygulama:	Water Quality