



Filter holder, for vacuum filtration, glass

Ürün Numarası:	234000
TRY Fiyat:	Bize Ulaşın
Mevcut	

Beaker-style funnel (300 mL) with 47 mm membrane support for filtration of aqueous samples. Fritted borosilicate glass membrane support for even filtration across surface of membrane. Funnel attaches evenly to base via aluminum clamp. Filter flask sold separately.

Teknik özellikler

Çap:	47 mm
Hacim / Paket boyutu:	1
Kutudakiler:	Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8).

Kutudakiler

Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8).